FitWidget

November 9, 2016

1 FitWidget

FitWidget is a widget to fit curves (1D data) with interactive configuration options, to set constraints, adjust initial estimate parameters...

1.1 Creating a FitWidget

First load the data.

```
In [1]: #%pylab inline
%gui qt

from silx.io import spech5


specfile = spech5.SpecH5("31oct98.dat")
xdata = specfile["/22.1/measurement/TZ3"]
ydata = specfile["/22.1/measurement/If4"]
#plot(xdata, ydata)
```

Then create a FitWidget.

```
In [2]: # Uncomment these lines if not in a jupyter notebook
    #from silx.gui import qt
    #a = qt.QApplication([])

from silx.gui.fit import FitWidget

fw = FitWidget()
  fw.setData(x=xdata, y=ydata)
  fw.show()
```



The selection of

fit theories and background theories can be done through the interface. Additional configuration parameters can be set in a dialog, by clicking the configure button, to alter the behavior of the estimation function (peak search parameters) or to set global constraints.

When the configuration is done, click the Estimate button. Now you may change individual constraints or adjust initial estimated parameters.


You can also add peaks by selecting *Add* in the dropdown list in the *Constraints* column of any parameter, or reduce the number of peaks by selecting *Ignore*.

When you are happy with the estimated parameters and the constraints, you can click the "Fit" button. At the end of the fit process, you can again adjust the constraints and estimated parameters, and fit again. Only click "Estimate" if you want to reset the estimation and all constraints (this will overwrite all adjustements you have done).

1.2 Open the FitWidget through a PlotWindow

A PlotWindow can be setup to show a fit icon. Clicking the icon will open a FitWidget with the active curve loaded as data.

```
In [3]: # Uncomment these lines if not in a jupyter notebook
     #from silx.gui import qt
     #a = qt.QApplication([])
```



FitConfig

```
from silx.gui.plot import PlotWindow

pw = PlotWindow(fit=True, control=True)
pw.addCurve(x=xdata, y=ydata)
pw.show()
```

All that is left to do is to select the curve and click the *Fit* icon. When the fitting is complete, the fitted function is plotted.

1.3 Exercice

Write a cubic polynomial function $y = ax^3 + bx^2 + cx + d$ and its corresponding estimation function (use a = 1, b = 1, c = 1, d = 1 for initial estimated parameters).

Generate synthetic data.

Create a FitWidget and add a cubic polynomial function to the dropdown list. Test it on the synthetic data.

1.3.1 Polynomial function

Tips: - Read the documentation for the module silx.math.fit.fittheory to use the correct signature for the polynomial and for the estimation functions. - Read the documentation for the module silx.math.fit.leastsq to use the correct format for constraints. Disable all constraints (set them to FREE)

Links: - http://pythonhosted.org/silx/modules/math/fit/fittheory.html#silx.math.fit.fittheory.FitTheory.fur - http://pythonhosted.org/silx/modules/math/fit/fittheory.html#silx.math.fit.fittheory.FitTheory.estimate

- http://pythonhosted.org/silx/modules/math/fit/leastsq.html#silx.math.fit.leastsq

1.3.2 Synthetic data

Tip: use the cubic_poly function

```
In [5]: import numpy
    x = numpy.linspace(0, 100, 250)
    a, b, c, d = 0.02, -2.51, 76.76, 329.14
    y = cubic_poly(x, a, b, c, d)
```

1.3.3 FitWidget with custom function

Tips:

you need to define a customized FitManager to initialize a FitWidget with custom functions

Doc:

http://www.silx.org/doc/silx/dev/modules/math/fit/fitmanager.html#silx.math.fit.fitmanager.FitManager.