

# Hdf5Widget

November 14, 2016

```
In [ ]: %gui qt
```

## 1 Required knowledges

- numpy (basic knowledge)
- Qt (basic knowledge)
- h5py (optionally)

## 2 Useful silx resources

- Getting start with the Hdf5 ([http://pythonhosted.org/silx/modules/gui/hdf5/getting\\_started.html](http://pythonhosted.org/silx/modules/gui/hdf5/getting_started.html))

## 3 Exercises summary

1. Features provided by Hdf5TreeView
  - Learn how to create an Hdf5TreeView
2. Create a HDF5 viewer
  - Learn how to use a dataset displayed by the tree
3. Create a tool to aggregate dataset and to create a diffraction mask
  - Use multi-selection node of the tree
  - Use features of the mask
4. Create a phase contrast viewer
  - Learn how to add context menu to the tree nodes

## 4 Features provided by Hdf5TreeView

```
from silx.gui import qt
from silx.gui import hdf5
app = qt.QApplication([])
tree = hdf5.Hdf5TreeView()
```

Name	Type	Shape	Value	Description	Node
nexus-20110325.h5					File
entry1					Group
entry2					Group
FE					Group
MM					Group
fe_v_gap	float64	141 × 1 ...			Dataset
end_time	string	1	['2011-03-25T13:50:45.025544']		Dataset
entry_identifier	string	1	['4']		Dataset
measurement					Group
fe_v_gap	float64	141 × 1 ...			Dataset
locum_l1	float64	141 × 1 ...			Dataset

HDF5 Tree

```
tree.setVisible(True)
app.exec_()
```

## 4.1 Exercise 0

1. Execute this script
2. Drag and drop an HDF5 file and play with it

## 4.2 Exercise 1

You can use `exercices/ex1_display.py` as skeleton

1. Create an application to load HDF5 file provided on the command line
  - Use [getting started with HDF5 widgets](#)

```
for filename in filenames:
    #
    # TODO: Load each filename into the model tree
    #
    print("Load %s" % filename)
```

## 4.3 Solution


```
In [ ]: !./working_examples/display.py nexus-20110325.h5
```

# 5 Create an HDF5 viewer

This exercise you how to use the `Hdf5TreeView` to browse and display datasets. We provide a `DataViewer` widget to help you to display the data.

## 5.1 DataViewer class

We provide a `DataViewer` widget, to display data using `Silx` plots.



HDF5 viewer

```
In [ ]: from silx.gui import plot
```

```
class DataViewer(qt.QStackedWidget):
    """Widget to display any kind of data"""

    def __init__(self, parent=None):
        """Constructor"""
        super(DataViewer, self).__init__(parent)

        self.__plot1d = plot.Plot1D()
        self.__plot2d = plot.Plot2D()
        self.__text = qt.QLabel()
        self.__text.setAlignment(qt.Qt.AlignCenter)

        self.__index1d = self.addWidget(self.__plot1d)
        self.__index2d = self.addWidget(self.__plot2d)
        self.__indexText = self.addWidget(self.__text)
        self.setCurrentIndex(self.__indexText)

    def showAsString(self, data):
        """Display a data using text"""
        self.__text.setText(str(data))
        self.setCurrentIndex(self.__indexText)

    def show1d(self, data):
        """Display a data using silx Plot1D"""
```

```

self.__plot1d.clear()
self.__plot1d.addCurve(legend="data", x=range(len(data)), y=data)
self.setCurrentIndex(self.__index1d)

def show2d(self, data):
    """Display a data using silx Plot2D"""
    self.__plot2d.clear()
    self.__plot2d.addImage(legend="data", data=data)
    self.setCurrentIndex(self.__index2d)

def show(self, data):
    """Display a data using the widget which fit the best"""
    isAtomic = len(data.shape) == 0
    isCurve = len(data.shape) == 1 \
        and numpy.issubdtype(data.dtype, numpy.number)
    isImage = len(data.shape) == 2 \
        and numpy.issubdtype(data.dtype, numpy.number)
    if isAtomic:
        self.showAsString(data)
    elif isCurve:
        self.show1d(data)
    elif isImage:
        self.show2d(data)
    else:
        self.showAsString(data)

```

Here is an example of use.

```

In [ ]: import numpy
        viewer = DataViewer()
        viewer.setVisible(True)

In [ ]: # To display an image
        viewer.show(numpy.random.rand(100, 100))

In [ ]: # or a curve
        viewer.show(numpy.random.rand(100))

In [ ]: # or a value
        viewer.show(numpy.random.rand(1)[0])

```

## 5.2 Viewer class

We also provide a Viewer class. This class display together an Hdf5TreeView and a DataViewer.

```

In [ ]: window = qt.QSplitter()
        tree = hdf5.Hdf5TreeView(window)
        viewer = DataViewer(window)

```

```

window.addWidget(tree)
window.addWidget(viewer)
window.setStretchFactor(1, 1)
window.setVisible(True)

```

### 5.3 Exercise 2

You can use `exercises/ex2_viewer.py` as skeleton


1. Connect the tree to the viewer together
  - Use [getting started with HDF5 widgets](#)

### 5.4 Solution

```
In [ ]: !./working_examples/viewer.py nexus-20110325.h5
```

## 6 Create an aggregation from diffraction acquisition

This exercise show how to configure and use the `Hdf5TreeView` with multi-selection. It will be used to compute an aggregation on images. The use case is an aggregation of diffraction acquisitions in order to create a better mask.



HDF5 diffraction mask

### 6.1 Creating an aggregation

A sum of many images can be done like that with `numpy`. It is not the better way to have the best contrast for a diffraction mask, but is enough for this exercise.

```
In [2]: import numpy
        a = numpy.random.rand(5, 5)
        b = numpy.random.rand(5, 5)
        c = numpy.random.rand(5, 5)
        aggregate = numpy.sum([a, b, c], axis=0)
```

## 6.2 Exercise 3

You can use `exercises/ex3_diffraction_mask.py` as skeleton

1. Configure the tree as multi-selectable
  - Use [QAbstractItemView documentation](#)
2. Aggregate selected datasets on `onTreeActivated`
3. Show the result in the viewer
4. With the GUI, use the mask tool to create a mask from aggregated images

## 6.3 Solution

```
In [ ]: !./working_examples/diffraction_mask.py ID22_ma2909_Ti37Nb_450_72h_1.h5
```

# 7 Create an phase contrast viewer

This exercise show how to use the `Hdf5TreeView` context menu to a custom use. The use case is the phase contrast acquisition, in order to display better images from the raw data. To correct this images, we have to remove a background and apply a flat field. We can use the context menu to identify this dataset from an HDF5 file. The exercise provides few functions to help the computation.

## 7.1 Provided functions


The computation of corrected images is done using this equation using `raw`, `flatfield`, and `background` information.

$$corrected = \frac{raw - background}{flatfield - background}$$

```
In [ ]: def computeCorrectedImage(self, raw):
        if self.flatfield is None:
            raise RuntimeError("Flatfield is not defined")
        if self.background is None:
            raise RuntimeError("Background is not defined")

        raw = numpy.array(raw, dtype=numpy.float32)
        flatfield = numpy.array(self.flatfield.value, dtype=numpy.float32)
        background = self.background.value
        return (raw - background) / (flatfield - background)

        def setBackground(self, dataset):
```



HDF5 phase contrast viewer

```
self.background = dataset
```

```
def setFlatField(self, dataset):
    self.flatfield = dataset
```

## 7.2 Exercise 4

You can use `exercises/ex4_phase_contrast.py` as skeleton

1. Register a callback function for the context menu of the tree
  - Use [getting started with HDF5 widgets](#)
2. Create action to the menu to use the hovered dataset as background of flatfield
  - Use [getting started with HDF5 widgets](#)
3. Try to compute the corrected image when an image is selected in the tree and show it in the viewer

## 7.3 Solution

```
In [ ]: !./working_examples/phase_contrast.py ID16B_diatomee.h5
```