

io

November 10, 2016

1 silx IO modules

silx currently mainly provides modules to read SPEC data files. In a future release, the number of supported data formats will increase significantly by adding an optional dependency to the *FabIO* library.

A SPEC file can be converted into a HDF5 file.

silx has modules and functions to save configuration data stored in Python dictionaries to various formats (HDF5, INI, JSON).

It also provides a `save1D` function to save curves (1D numpy arrays) to simple file formats.

1.1 SPEC files

There are currently 3 modules handling SPEC file reading:

- `silx.io.specfile`: Low-level wrapping of the legacy C library
- `silx.io.specfilewrapper`: Second layer of wrapping to offer compatibility with the legacy PyMca wrapper
- `silx.io.spech5`: High-level module exposing a Spec file like a *h5py* file.

The use of `specfilewrapper` should be strictly limited to the case of existing programs using the PyMca wrapper, to switch to silx with minimal work: replacing the `import` statement `from PyMca5.PyMcaIO import specfilewrapper` with `from silx.io import specfilewrapper`.

New programs should use `specfile` or (preferably) `spech5`.

We encourage using of `spech5`, because it provides the same interface that *h5py* uses for reading HDF5 files. The stated goal of *silx* regarding IO is to provide a unified API to as many data formats as possible, to simplify training.

A fourth module related to SPEC file deals with conversion of SPEC files into HDF5 files: `silx.io.spectoh5`

1.1.1 spech5

Documentation: <http://www.silx.org/doc/silx/dev/modules/io/spech5.html>

Exposed structure

```
/
  1.1/
    title = "...
    start_time = "...
    instrument/
      specfile/
        file_header = "...
        scan_header = "...
      positioners/
        motor_name = value
        ...
      mca_0/
        data = ...
        calibration = ...
        channels = ...
        preset_time = ...
        elapsed_time = ...
        live_time = ...

      mca_1/
        ...
      ...
    measurement/
      colname0 = ...
      colname1 = ...
      ...
      mca_0/
        data -> /1.1/instrument/mca_0/data
        info -> /1.1/instrument/mca_0/
        ...
  2.1/
    ...
```

Opening a Spec file and listing all scans Like a *h5py* file, a SPEC file is composed of a tree structure of groups (~folders) and datasets (~files).

The root group contains one subgroup per scan.

```
In [1]: from silx.io import spech5
```

```
sfh5 = spech5.SpecH5("31oct98.dat")
```

```
#print a list of scans
print(sfh5.keys())
```

```
['19.1', '20.1', '21.1', '22.1', '23.1', '24.1', '25.1', '26.1', '27.1', '28.1', '2
```

Accessing detector measurements Measurements are 1D datasets (data columns in a scan in the original file). They are located in a sub-group of your scan named “measurement”.

You can iterate over all keys in a SpecH5Group. Each key is the name of a subgroup or a dataset.

```
In [2]: # print all detector labels in all scans
        for scan_key in sfh5:
            scan_group = sfh5[scan_key]
            print("####scan number " + scan_key + "####")
            if scan_key.startswith("6"):
                print(scan_group["measurement"].keys())
```

```
####scan number 19.1####
####scan number 20.1####
####scan number 21.1####
####scan number 22.1####
####scan number 23.1####
####scan number 24.1####
####scan number 25.1####
####scan number 26.1####
####scan number 27.1####
####scan number 28.1####
####scan number 29.1####
####scan number 30.1####
####scan number 31.1####
####scan number 32.1####
####scan number 33.1####
####scan number 34.1####
####scan number 35.1####
####scan number 36.1####
####scan number 37.1####
####scan number 38.1####
####scan number 39.1####
####scan number 40.1####
####scan number 41.1####
####scan number 42.1####
####scan number 43.1####
####scan number 44.1####
####scan number 45.1####
####scan number 46.1####
####scan number 47.1####
####scan number 48.1####
####scan number 49.1####
####scan number 50.1####
####scan number 51.1####
####scan number 52.1####
####scan number 53.1####
####scan number 54.1####
```

```

####scan number 55.1####
####scan number 56.1####
####scan number 57.1####
####scan number 58.1####
####scan number 59.1####
####scan number 60.1####
['TX3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 61.1####
['TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 62.1####
['TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 63.1####
['TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 64.1####
['TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 65.1####
['TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 66.1####
['TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 67.1####
['TX3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 68.1####
['TX3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 69.1####
['TX3', 'TZ3', 'Epoch', 'Seconds', 'If2', 'If3', 'If5', 'If6', 'If7', 'If8', 'IO', 'It', 'V']
####scan number 678.1####
['Bal', 'Blo', 'Bli']
####scan number 679.1####
['Bal', 'Blo', 'Bli']
####scan number 680.1####
['Bal', 'Blo', 'Bli']
####scan number 681.1####
['Bal', 'Blo', 'Bli']
####scan number 682.1####
['Bal', 'Blo', 'Bli']
####scan number 684.1####
['Bal', 'Blo', 'Bli']
####scan number 685.1####
['Bal', 'Blo', 'Bli']
####scan number 687.1####
['Bal', 'Blo', 'Bli']
####scan number 688.1####
['Bal', 'Blo', 'Bli']
####scan number 700.1####
####scan number 688.2####
['Bal', 'Blo', 'Bli']
####scan number 688.3####
['Bal', 'Blo', 'Bli']

```

Datasets are objects that are similar to numpy arrays (or h5py datasets). You can access the individual values using indices, or extract subarrays using slicing.

```
In [3]: # get two data columns
        xdata = sfh5["/22.1/measurement/TZ3"]
        ydata = sfh5["/22.1/measurement/If4"]


        print(xdata)
        print(sfh5["/22.1/measurement/TZ3"][5:8])
```

29.29999924	29.30397606	29.30795288	29.31202888	29.3160038	
29.31998062	29.32395554	29.32803154	29.33200836	29.33598328	
29.3399601	29.3440361	29.34801292	29.35198784	29.35596466	
29.36004066	29.36401558	29.3679924	29.37196732	29.37604332	
29.38002014	29.38399506	29.38797188	29.39204788	29.3960247	
29.39999962	29.40397644	29.40795135	29.41202736	29.41600418	
29.4199791	29.42395592	29.42803192	29.43200874	29.43598366	
29.43996048	29.44403648	29.4480114	29.45198822	29.45596313	
29.46003914	29.46401596	29.46799278	29.4719677	29.4760437	
29.48002052	29.48399544	29.48797226	29.49204826	29.49602318	29.5
29.50397682	29.50795174	29.51202774	29.51600456	29.51997948	
29.5239563	29.5280323	29.53200722	29.53598404	29.53996086	
29.54403687	29.54801178	29.5519886	29.55596352	29.56003952	
29.56401634	29.56799126	29.57196808	29.57604408	29.5800209	
29.58399582	29.58797264	29.59204865	29.59602356	29.60000038	
29.6039753	29.60795212	29.61202812	29.61600494	29.61997986	
29.62395668	29.62803268	29.6320076	29.63598442	29.63995934	
29.64403534	29.64801216	29.65198708	29.6559639	29.6600399	
29.66401672	29.66799164	29.67196846	29.67604446	29.68001938	
29.6839962	29.68797112	29.69204712	29.69602394	29.70000076	
29.70397568	29.7079525	29.7120285	29.71600342	29.71998024	
29.72395515	29.72803116	29.73200798	29.7359848	29.73995972	
29.74403572	29.74801254	29.75198746	29.75596428	29.76004028	
29.7640152	29.76799202	29.77196884	29.77604485	29.78001976	
29.78399658	29.7879715	29.7920475	29.79602432	29.79999924	
29.80397606	29.80795288	29.81202888	29.8160038	29.81998062	
29.82395554	29.82803154	29.83200836	29.83598328	29.8399601	
29.8440361	29.84801292	29.85198784	29.85596466	29.86004066	
29.86401558	29.8679924	29.87196732	29.87604332	29.88002014	
29.88399506	29.88797188	29.89204788	29.8960247	29.89999962]	
[29.31998062	29.32395554	29.32803154]			

```
In [4]: %pylab inline
        plot(xdata, ydata)
```

Populating the interactive namespace from numpy and matplotlib

```
Out[4]: [<matplotlib.lines.Line2D at 0x7f9b5a2f3a20>]
```




Accessing MCA data mca data can be found in /94.1/measurement/mca_0/data or /94.1/instrument/mca_0/data. This is a 2D array with all spectra recorded for one analyser during a scan. If there are more than one analyser, you will find more groups (mca_1, mca_2...)

Scan metadata (headers) can be found in mca_0/info.

```
In [5]: sfh5_2 = spech5.SpecH5("oleg.dat")
        scan94_mca0 = sfh5_2["94.1/measurement/mca_0"]
        print(scan94_mca0.keys())
        mca_data = scan94_mca0["data"]
        print(mca_data.shape)
        for i, mca_spectrum in enumerate(mca_data[10:15]):
            plot(mca_spectrum[475:550], label="spectrum %d" % (10 + i))
        legend()

['data', 'info']
(21, 1024)
```

```
Out[5]: <matplotlib.legend.Legend at 0x7f9b5c349e10>
```



```
In [6]: mca_info = scan94_mca0["info"]
        print(mca_info.keys())
        for datasatet_name in mca_info:
            print("dataset %s: %s" % (datasatet_name, mca_info[datasatet_name]))

['data', 'calibration', 'channels']
dataset data: [[ 0.  0.  0. ...,  0.  0.  0.]
 [ 0.  0.  0. ...,  0.  0.  0.]
 [ 0.  0.  0. ...,  0.  0.  0.]
 ...,
 [ 0.  0.  0. ...,  0.  0.  0.]
 [ 0.  0.  0. ...,  0.  0.  0.]
 [ 0.  0.  0. ...,  0.  0.  0.]]
dataset calibration: [ 0.  1.  0.]
dataset channels: [  0   1   2 ..., 1021 1022 1023]
```

1.2 save1D

The formats supported by save1D are:

- .csv
- .npy (format used by `numpy.save` / `numpy.load`)
- .txt (CSV with customized header, footer and record separator)
- .dat (simplified SPEC file, one scan with two columns per curve)

For example, you can save the data loaded previously as txt:

```
In [7]: from silx.io.utils import save1D

        hdr = "This is data extracted from a Spec file"
        ftr = "My custom footer"

        save1D("If4_TZ3.txt", x=xdata, y=ydata, xlabel="TZ3", ylabels=["If4"], csv=
```

y can be a 2D array in case you want to save multiple curves to a single file. That's why ylabels is formatted as a list of labels.

1.3 config dict

Python dictionaries are very convenient to handle configuration data for a program. The silx.io.dictdump and silx.io.configdict modules address the issue of dumping dictionaries to files and loading them back.

```
In [9]: from silx.io.dictdump import load, dump
        import numpy

        mydic0 = {
            'simple_types': {
                'float': 1.0,
                'int': 1,
                'string': 'Hello World',
            },
            'containers': {
                'list': [-1, 'string', 3.0, False],
                'array': numpy.array([1.0, 2.0, 3.0]),
                'dict': {
                    'key1': 'Hello World',
                    'key2': 2.0,
                }
            }
        }

        # dump(mydic0, "mydic.json")
        # dump(mydic0, "mydic.h5")
        dump(mydic0, "mydic.ini")

In [10]: # mydic1 = load("mydic.json")
        # mydic1 = load("mydic.h5")
        mydic1 = load("mydic.ini")

        print(mydic1)
```

```
ConfigDict([('simple_types', OrderedDict([('string', 'Hello World'), ('int', 1), ('float', 1.0), ('key2', 2.0), ('key1', 'Hello World')])])])
```


The JSON format cannot serialize numpy arrays. Except for this drawback, it is a very common and convenient ascii format to handle simple python data.

The HDF5 format handles numpy arrays very well, but data types are not preserved in heterogeneous lists (lists containing a mix of strings, integers...), all elements are cast into byte strings for these lists.

The INI format handles all data types well in our example. Its drawback is that the INI format is not standardized, and the INI files produced by `silx.configdict`, that are designed to handle nested configuration dictionaries, are not commonly read by other libraries than `silx`.

2 spectoh5

The `silx.io.spectoh5` module provides 2 functions to convert SPEC files to HDF5 files: - `convert()`: simple function to convert a file - `write_spec_to_h5()`: function with more customization options, such as specifying a target path/group inside the output HDF5 file where scans will be written.

The documentation for these functions is: - <http://www.silx.org/doc/silx/dev/modules/io/spectoh5.html#>
- http://www.silx.org/doc/silx/dev/modules/io/spectoh5.html#silx.io.spectoh5.write_spec_to_h5

2.1 Exercise

1. Convert a SPEC file into an HDF5 file, and list it's content with `silx.io.utils.h5ls`
2. Write two different SPEC files into a single HDF5 file, into 2 separate target groups.
3. Redo question 1. or 2., with auto-chunking and gzip compression

Tips: - documentation for `h5ls`: <http://www.silx.org/doc/silx/dev/modules/io/utils.html#silx.io.utils.h5ls>
- additional documentation for question 3.: + http://docs.h5py.org/en/latest/high/group.html#Group.create_dataset
+ <http://docs.h5py.org/en/latest/high/dataset.html#dataset-compression>

```
In [ ]: # Question 1.
```

```
from silx.io.spectoh5 import convert
from silx.io.utils import h5ls

convert("oleg.dat", "oleg.h5", mode="w")
print(h5ls("oleg.h5"))
```

```
In [ ]: # Question 2.
```

```
from silx.io.spectoh5 import write_spec_to_h5
from silx.io.utils import h5ls

# specifying a path causes the file to be open,
# then closed again after writing is done
write_spec_to_h5("31oct98.dat", "concat.h5",
                 h5path="/31oct98", mode="w")
write_spec_to_h5("oleg.dat", "concat.h5",
                 h5path="/oleg", mode="a")
```

```
# alternative way: use a file handle rather than a file path
# (you are responsible for opening/closing the file)
import h5py

output_file = h5py.File("concat.h5", mode="a")
# write_spec_to_h5("31oct98.dat", output_file,
#
h5path="/31oct98", overwrite_data=True)
write_spec_to_h5("oleg.dat", output_file,
                h5path="/oleg_2")
output_file.close()
```

```
In [ ]: # Question 3.
from silx.io.spectoh5 import convert
from silx.io.utils import h5ls

create_ds_args = {"chunks": True,
                  "compression": "gzip"}

convert("oleg.dat", "oleg.h5", mode="w",
        create_dataset_args=create_ds_args)
print(h5ls("oleg.h5"))
```

```
In [ ]:
```