

Input/output


Silx IO

- Provides function to open file data
 - Accessible as HDF5-like objects
 - In read-only
 - Support HDF5 files, Spec files, EDF files (plus format supported by *FabIO*)
- In symbiosis with our widgets
 - HDF5 tree widget
 - DataViewer
- Also contains
 - Maintained *specfilewrapper* class
 - Maintained TIFF and EdfFile reader and writer
 - Spec to HDF5 converter
 - Dictionary dump

HDF5 introduction

HDF5 (for Hierarchical Data Format) is a file format to structure and store data.

- Standard exchange format for heterogeneous data
- Hierarchical collection of data (directory and file, UNIX-like path)
 - **File:** the root of the container
 - **Group:** a grouping structure containing groups or datasets
 - **Dataset:** a multidimensional array of data elements
 - And other features (links, attributes, datatypes)



HDF5 example

Here is an example of file generated by pyFAI.


Name	Type	Shape	Value	Description	Node
test.h5					File
diff_map_0000					Group
data					Group
map	float32	29 × 78 × 100 ...			Dataset
program_name	string		pyFAI		Dataset
pyFAI					Group
PONIfile	string		None		Dataset
date	string		2016-06-09T14:03:18+02:00		Dataset
detector	string		Pilatus1M		Dataset
dim0	int64		29		Dataset
dim1	int64		78		Dataset
dim2	int64		100		Dataset
dist	float64		1.63410088403		Dataset
inputfiles	string	2230	['/home/valls/workspace/data/diff_ma...		Dataset
pixel1	float64		5.197559e-05		Dataset
pixel2	float64		5.158538e-05		Dataset
poni1	float64		0.0312319583221		Dataset
poni2	float64		0.0141040323735		Dataset
program	string	1	['diff_map']		Dataset
rot1	float64		-0.0102999203254		Dataset
rot2	float64		0.00868154231655		Dataset

NeXus introduction

NeXus is a common data format for neutron, x-ray and muon science based on HDF5.

Set of rules to structure a HDF5 file for interoperability of the data.


We try to follow NeXus convention to display data or to expose data.



Specfile with silx

Silx provides access to spec files using an HDF5-like mapping. It uses a subset of the HDF5 model, based on NeXus.


HDF5-like mapping



EDF with silx

Silx provides access to EDF files (and other format supported by *FabIO*) using an HDF5-like mapping. It is a subset of the HDF5 model, based on NeXus.

HDF5-like mapping



Silx IO API

Based on HDF5 and *h5py* module API.

Open a file

```
import silx.io

# open
obj = silx.io.open(filename)
# do your stuff here
obj.close()

# or using context manager
with silx.io.open(filename) as obj:
    # do your stuff here
    # the close is called for you at the end of the with
```


Silx IO API

Based on HDF5 and *h5py* module API.

Common properties

```
obj.name    # the path name
obj.parent  # the direct container of the object
obj.file    # the file container of the object

# test object type
if silx.io.is_file(obj):
    print("this is a root file")

    # path of the file from the file system
    obj.filename

if silx.io.is_group(obj):
    # BTW a file is a group
    print("this is a group")

if silx.io.is_dataset(obj):
    print("this is a dataset")
```

Silx IO API

Based on HDF5 and *h5py* module API.

Node traversal

```
if silx.io.is_group(obj):
    # it contains child

    # number of child
    len(obj)

    # iterator on child names
    obj.keys()

    # access to a child
    child = obj["child_name"]

    # access to a child using a path
    child = obj["path/to/a/child"]

    # the path can be absolute
    child = obj["/absolute/path/to/a/child"]
```

Silx IO API

The content of a dataset is a *numpy* data.

Data access

```
if silx.io.is_dataset(obj):
    # it contains data

    # a dataset provides information to the data
    obj.shape    # multidimensional shape
    obj.size     # amount of items
    obj.dtype    # type of the array

    # copy the full data as numpy array
    data = obj[...]

    # or a part of it (using numpy selector)
    data = obj[1:2, ::3, 7]

    # special case to access to the value of a scalar
    # i.e. a single integer, a single string...
    data = obj[()]
```

Silx IO API

Specfile example

```
import silx.io

h5like = silx.io.open('data/oleg.dat')

# print available scans
print(h5like['/'].keys())

# print available measurements from the scan 94.1
print(h5like['/94.1/measurement'].keys())

# get data from measurement
time = h5like['/94.1/measurement/Epoch']
bpm = h5like['/94.1/measurement/bpmi']
mca = h5like['/94.1/measurement/mca_0/data']
```

Silx IO API

EDF example

```
import silx.io.utils

h5like = silx.io.open("data/medipix.edf")


# here is the data as a cube using numpy array
# it's a cube of images * number of frames
data = h5like["/scan_0/instrument/detector_0/data"]
# here is the first image
data[0]

# groups containing datasets of motors, counters
# and others metadata from the EDF header
motors   = h5like["/scan_0/instrument/positioners"]
counters = h5like["/scan_0/instrument/measurement"]
others   = h5like["/scan_0/instrument/detector_0/others"]

# reach a monitor named 'mon'
# it's a vector of values * number of frames
monitor = counters["mon"]
# here is the monitor value at the first frame
monitor[0]
```

HDF5 tree

Tree view which allow to browse HDF5 file content.



The screenshot shows a window titled 'HDF5 tree' with a table of file contents. The table has columns: Name, Type, Shape, Value, Description, and Node. The tree structure is as follows:

Name	Type	Shape	Value	Description	Node
nexus-20110325.h5	File				File
entry1	Group				Group
entry2	Group				Group
FE	Group				Group
MM	Group				Group
fe_v_gap	Dataset	float64 141 × 1 ...			Dataset
end_time	Dataset	string 1	['2011-03-25T13:50:45.025544']		Dataset
entry_identifier	Dataset	string 1	['4']		Dataset
measurement	Group				Group
fe_v_gap	Dataset	float64 141 × 1 ...			Dataset
locum_l1	Dataset	float64 141 × 1 ...			Dataset

Provides:

- Mouse click event
- Customable content menu
- Customable columns

HDF5 tree

Example


Tree view which allow to browse HDF5 file content.

```
import silx.gui.hdf5
tree = silx.gui.hdf5.Hdf5TreeView()
model = tree.findHdf5TreeModel()

# Insert a filename
model.insertFile("data/test.h5")

# Insert an HDF5 node
h5 = silx.io.open("data/test.h5")
model.insertH5pyObject(h5)
```

DataViewer



View of the
selected data

Selection of axes
and slicing

Selection of the
viewer mode

DataViewer

Example

```
import silx.gui.data.DataViewerFrame
viewer = silx.gui.data.DataViewerFrame.DataViewerFrame()
viewer.setVisible(True)

# Let's display a random cube
import numpy
data = numpy.random.rand(100, 100, 100)
viewer.setData(data)

# Let's display an HDF5 dataset
import silx.io
h5like = silx.io.open("data/ID16B_diatomee.h5")
dataset = h5like["/data/0299"]
viewer.setData(dataset)
```

Exercises

This exercises are based on phase contrast acquisition data.

You can find it as notebook, or as Python files.

- **Exercise 1:**
 - Browse the data file
- **Exercise 2:**
 - Compute the correction
- **Exercise 3:**
 - Create a viewer
- **Exercise 4:**
 - Custom the viewer to display corrected images